Osteo-Comfort

Gebrauchsanweisung User manual

◆ Inhalt

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Content

- 8 Product Description
- 8 Indication
- 8 Contraindication
- **8** Function
- 8 Reimbursement Code
- 9 Area of Application
- 9 REF/Sizes
- **10** Safety
- **11** Application
- 13 Comfort/Hygiene
- 13 Storage/Disposal
- **13** Material
- **13** Product Care/Cleaning/Maintenance

Product description

Osteo-Comfort is a highly functional back orthosis, which consists of a spine reclinator and a circular lumbar orthosis. It dynamically straightens the body and thus promotes muscle activation.

◆ Indication

- Stable vertebrate fractures
- Spinal compression fracture (stable fracture)
- Lumbar instability
- Impression fracture
- Spinal tumor

Contraindication

Non-compliant patients

Function

- Straightening of the shoulder girdle
- Stabilisation of the spine
- Support of the dorsal and abdominal muscles
- Combines a spine reclinator with a circular lumbar orthosis
- It is possible to don the orthosis without twisting the upper body
- Closing requires only little strength thanks to the six-fold pulley system

Reimbursement Code

DE: HMV-Nr.: 23.15.04.2004 AT: Tarif-Nr.: 31013 026

Application Area

• TH8 - L5*

*In general, the orthosis should only be worn by patients with stable vertebral fractures. The application area given above is an application recommendation arising from the biomechanical function of the orthosis. An individual evaluation by a doctor will always be necessary to determine whether the orthosis is suited for the patient in question.

REF/Sizes

Osteo-Comfort Orthosis, long			
REF	Size	Girth	Body Height
22.610.020.00	XS	55 - 70 cm	over 165 cm
22.610.021.00	S	65 - 80 cm	over 165 cm
22.610.022.00	М	75 - 90 cm	over 165 cm
22.610.023.00	L	85 - 100 cm	over 165 cm
22.610.024.00	XL	95 - 110 cm	over 165 cm
22.610.025.00	2XL	105 - 120 cm	over 165 cm
22.610.026.00	3XL	115 - 130 cm	over 165 cm

Osteo-Comfort Orthosis, short			
REF	Size	Girth	Body Height
22.610.010.00	XS	55 - 70 cm	under 165 cm
22.610.011.00	S	65 - 80 cm	under 165 cm
22.610.012.00	M	75 - 90 cm	under 165 cm
22.610.013.00	L	85 - 100 cm	under 165 cm
22.610.014.00	XL	95 - 110 cm	under 165 cm
22.610.015.00	2XL	105 -120 cm	under 165 cm
22.610.016.00	3XL	115 - 130 cm	under 165 cm

Safety

- Before this device is used, it is essential that the orthosis is individually adapted
 by an experienced orthopedic technician and the patient is trained in the orthosis' use according to the doctor's instructions.
- Incorrect selection, adaptation, and/or failure to monitor the correct fit of this product may cause health damage.
- Never wear the orthosis on broken skin.
- Check orthosis for damage (such as loose connections, fasteners, belts, etc.) each time it is put on.
- Should any performance malfunctions (like an apparent reduction or lack of effect) be suspected, immediately call on your doctor or orthopedic technician.
- The orthosis must be checked for proper function, correct fit and possible damages by a doctor or orthopedic technician at regular intervals. The timing of these intervals has to be defined for each patient individually.
- Should the product have been subject to disproportionally strong forces, it must be checked for possible damages by an orthopedic technician before further use.
- Harmful pressure caused by the orthosis may cause skin injuries or disrupt the blood flow and must be avoided

Application

1. Scope of Delivery:

- Straps (including shoulder pad)
- Back pad
- Reclinator
- Basic component for lumbar orthosis
- Enclosed hook and loop strap

2. Attaching the Shoulder Straps
Open the zipper on the pad. Pull the recliner slightly out of the pocket.
Slide the strap ends out of the reclinator (towards the drillhole).
Pull the separated straps out of the pocket. (a) Completely remove the seperated reclinator.

3. Attaching the Reclinator to the Basic Component of the Lumbar Orthosis

Connect the reclinator with the circular orthosis using the quick-change fasteners (components can be snapped together). To seperate them, undo the fasteners with a screwdriver.

4. The degree to which the upper body is straightened can be adjusted by repositioning the strap of the lumbar orthosis. Slight straightening effect = front buttonhole / strong straightening effect = back buttonhole

5. To make the orthosis easier to use for your patient, place the enclosed piece of hook and loop fastener onto the center of the reclinator, about 5 cm above the bottom end of the zipper. Then attach the two ends to the inside of the left and the right side of the basic lumbar component. The basic component will now be held in place when the orthosis is donned, and putting the orthosis on will be a lot easier.

 Open the Velcro fastening of the circular orthosis. Put your arms through the shoulder straps, so the orthosis sits on your back like a backpack.

The length of the shoulder straps can be adjusted in two directions. Adjust the length of the strap which is going upwards first. If it is necessary to shorten the strap, lead the strap end towards the shoulder pad, and attach it on the strap, using the hook and loop fastener.

 Pull both ends of the circular orthosis downward/forward. Now close the Velcro fastening of the circular orthosis in front of the belly.

8. Using the handle, draw the pulley forwards, until you can feel the compression. Attach the handle to the circular orthosis with the velcro patch.

9. Pull the strap ends down/towards the back until you can feel a straightening effect.

Lead the ends of the strap which has been pulled downwards back up, and attach them to the strap using the hook and loop fastener.

Comfort/Hygiene

Due to the required individual adaptation (product-dependent forming and/ or machining) of the orthosis it is not suitable for reuse by other persons!

Storage/Disposal

Product does not contain any environmentally hazardous substances. It can be disposed of with household refuse.

Material

Textile components (latex free):

- 50% Polyester
- 48% Polyurethane
- 2% Lycra

Reclinator: Spring steel

◆ Product Care/ Cleaning/ Maintenance

- Hand wash at 30°C
- Do not iron
- Do not tumble dry
- Do not bleach
- Do not dry clean

◆ Notizen / Notes	

◆ Notizen / Notes	

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Wilhelm Julius Teufel GmbH Robert-Bosch-Straße 15 73117 Wangen Deutschland